

THE ONYFIX® NAIL CORRECTION SYSTEM: PAINLESS – SAFE – EFFECTIVE – INNOVATIVE

The Onyfix[®] Nail Correction System is a revolutionary new product for the treatment of ingrown and involuted nails. Your foot specialist will apply this customized system onto the affected nail, which then corrects the shape over time through natural nail growth. The system is completely painless to apply and to wear!



THE REVOLUTIONARY PAIN-FREE, NON-INVASIVE NAIL CORRECTION SYSTEM

QUESTIONS? PLEASE CONTACT US, WE'RE HAPPY TO HELP!

Make an appointment today with our Certified Onyfix® Specialist!

Onyfix[®] services are available at:

O Painless application process

- Non-invasive no needles or surgery
- O Provides immediate, lasting pain relief
- O Corrects the nail shape through natural nail growth
- O Evidence-based studies proving efficacy
- () No restrictions on activities (swimming, sauna, nail polish, etc.)
- Suitable for all, including children, high-risk clients and those with diabetes



neubourg healthcare Inc. | 1500 - 480 University Avenue Toronto, ON M5G 1V2 | Canada | www.onyfix.com



8000_02

GET HEALTHY, BEAUTIFUL NAILS AGAIN WITH ONYFIX®! HOW DOES ONYFIX® WORK?

NEV

The hardened Onyfix® strip ensures that the shape is maintained during your nail's natural

growth. Some nails might require more than one

cycle in order to properly correct the curvature.

After successful treatment, your foot specialist

will remove the Onyfix® from your nail easily and

PAINFUL, INGROWN **TOENAILS?**

Did you know that:

- O Ingrown toenails that keep growing untreated can grow deeper and deeper into the tissue?
- O Ingrown nails can cause pain with every single step and even change the way you walk?
- O Inflammation of the nail fold can cause serious complications?

Don't continue suffering with pain or risk serious complications from untreated ingrown nails. Contact your Certified Onyfix[®] Specialist today for your painless, customized Onyfix® application!

TREATMENT TIMELINE

require protective eyewear!

(1.)

2.)



Start

After 10 weeks

After 20 weeks

(3.)

(4.)

painlessly!

After 30 weeks



With Onyfix[®], almost all types of ingrown and involuted nails are treatable. Your certified Onyfix[®] specialist will first roughen

and clean the nail in preparation for the application.

The Onyfix[®] is like a soft putty and is applied to the

nail. A spatula is used to shape it into a strip and

is then used to harden the Onyfix® and bond it to

the nail. Don't worry, the light is safe and doesn't

place it in the right position. A special blue light